



CONTOUR

NON-INVASIVE  
FACE & BODY  
CONTOURING

3D HIGH INTENSITY  
FOCUSED ULTRASOUND





# MÖ CONTOUR

SCANDINAVIAN  
SKINCARE SYSTEMS



AIMING LIGHT



MAIN CONSOLE



HANDLE CONNECTOR



HANDPIECE



HANDLE SUPPORT



NON-INVASIVE  
ONE TREATMENT NEEDED  
TIGHTENS  
YOUTHFUL LOOK  
FAST OPERATION TIME  
PAINLESS

NO DOWN TIME  
ACHIEVES INSTANT RESULTS  
REJUVENATES  
NO SIDE EFFECTS  
LONG LASTING RESULTS  
SAFE





# THE TECHNOLOGY

MÖ CONTOUR IS THE LATEST PIECE OF EQUIPMENT IN THE MÖ COLLECTION RANGE.

USING NON-INVASIVE LIFTING AND TIGHTENING TECHNOLOGY, IT REJUVINATES AND REDUCES SIGNS OF AGING.

MÖ CONTOUR USES 3D HIGH INTENSITY FOCUSED ULTRASOUND TO ACTIVATE THE SKIN'S THERMAL COAGULATION PROCESS, THEREBY PROVIDING INSTANT AND LONG-LASTING RESULTS.

MÖ CONTOUR IS AN EXTREMELY EFFECTIVE TREATMENT, USING SAFE AND PAINLESS PROCEDURES.





# THE SPEC

TYPE	3D HIGH INTENSITY FOCUSED ULTRASOUND
FREQUENCY	4MHZ
POWER	200W
DEPTH (STANDARD)	DL 1.5MM, DL 3.0MM, DL 4.5MM
DEPTH (OPTIONAL)	DL 6.0MM, DL 8.0MM, DL 10.0MM, DL 13.0MM, DL 16.0MM
LENGTH	5 - 25MM
WIDTH	0 - 10MM
ENERGY	0.2 - 2.0J
PITCH	1.0MM - 5.0MM
CLEARANCE	1.0MM - 5.0MM
POWER SOURCE	AC 110V ~ 220V ~ 50Hz / 60Hz
MACHINE SIZE & NET WEIGHT	400 X 350 X 335MM & 8.5Kg
PACKAGE SIZE & GROSS WEIGHT	525 X 445 X 465MM & 12Kg
TROLLEY SIZE & NET WEIGHT	540 X 480 X 790MM & 10Kg
TROLLEY PACKAGE SIZE & GROSS WEIGHT	680 X 680 X 330MM & 15Kg





## WHO IS SUITABLE FOR MÖ CONTOURING?

THE PROGRAMME IS SUITABLE FOR CLIENTS AGED BETWEEN 30 AND 60, WHO WANT TO REDUCE PORES, TIGHTEN THE SKIN AND LOOK YOUNGER.

## WHEN WILL YOU SEE RESULTS AND WHAT ARE THE SIDE EFFECTS?

80% OF RESULTS ARE SEEN AFTER RECOVERY DEPENDING ON THE CLIENT'S SKIN PROPERTY AND CELL- REPAIR (SUCH AS DPHP). CLIENTS OVER 30 YEARS HAVE SIGNIFICANTLY LESS COLLAGEN AND WILL BE RECOMMENDED TO DRINK SPECIAL REPAIR OF COLLAGEN ACMETEA BEFORE TREATMENT. THIS WILL HELP REGENERATION OF COLLAGEN FIBERS. IF CLIENTS DO NOT HAVE ENOUGH COLLAGEN, DUE TO DEFICIENCY OF CELL REGENERATION IT CAN RESULT IN HOLLOW CHEEK.

## WHAT IS THE SUITABLE INTERVAL TIME?

USUALLY THE PROGRAMME RESULTS LAST UP TO 2-3 YEARS. FACE AND BODY CONTOURING CAN DAMAGE THE SUBCUTANEOUS, SO AN INTERVAL TIME OF 2 YEARS IS RECOMMENDED.

## WHAT IS THE POST CARE ADVICE?

AFTER THE TREATMENT, THE PATIENT CAN GO ABOUT DAILY LIFE AS THERE ARE NO SPECIAL CONSIDERATIONS. THE TREATMENT AREA MAY BE A LITTLE SWOLLEN, BUT THIS WILL DISAPPEAR AFTER A FEW HOURS.



# THE HANDPIECE

SLEEK & FUNCTIONAL DESIGN

HIGH FOCUS TARGETING

LATEST TECHNOLOGY

HIGH QUALITY


SHAPED FOR COMFORT

EASY TO CLEAN

1 YEAR GUARANTEE

SAFETY ASSURED

CUSTOMISABLE TO AREA







DL 1.5<sub>MM</sub>  
FULL FACE  
& NECK



DL 3.0<sub>MM</sub>  
FULL FACE  
& NECK



DL 4.5<sub>MM</sub>  
FULL FACE  
& NECK



DL 6.0<sub>MM</sub>  
BACK, CHEST  
& ABDOMEN

T  
R  
E  
A  
T  
M  
E  
N  
T  
  
G  
U  
I  
D  
E

DL 8.0<sub>MM</sub>

BACK, CHEST  
& ABDOMEN



DL 10<sub>MM</sub>

BACK, LEG  
& ABDOMEN



DL 13<sub>MM</sub>

LEGS &  
ABDOMEN



DL 16<sub>MM</sub>

LEGS &  
ABDOMEN

