



PLASMAFILL



The Next Generation of Filler

The MÖ PlasmaFill is an autologous method for medical rejuvenation using 'bio-filler', and is a natural alternative to hyaluronic acid fillers.

PlasmaFill is a holistic approach to promote a natural volume for skin rejuvenation, and the product is easily absorbed by the body. The process uses the waste blood product (platelet poor plasma- PPP) from the procedure 'Platelet Rich-Plasma' (PRP) - using your own blood's plasma.

This is an alternative to synthetic fillers for those who are allergic to HA products or are just not keen on the HA aesthetic look.

The MÖ PlasmaFill can be used for:

- Soft filler
- Deep filler
- Skin-boosting, 5 & 8 points.
- Volume
- Smoothing out deep lines & wrinkles
- Reducing fine lines (especially around the face, neck and body)
- Enhancing youthful appearance

This is a much safer alternative to dermal fillers and the results provided are typically even better than you get with dermal filler. One major benefit of using bio-filler as a cosmetic filler is its ability to help your skin heal and rejuvenate faster, without the risk of an allergic reaction or skin sensitivity. If PRP is accidentally injected into a vessel under your eye, the risk of blindness is essentially zero.

It is 100% safe since it uses your own blood's plasma and has no dangerous side effects.

What are the side effects of PlasmaFill?


Generally patients can expect some light bruising, redness and swelling around where the injections were administered, along with a little tenderness and pain there.

These should disappear within a few days and is nothing abnormal.

How long does it last?

The MÖ PlasmaFill treatment can last from 3 - 4 months when using soft filler with thin density and 8 - 10 months with thick filler after one session.

Results can last over a year when using the skin-boosting, 5 & 8 point method after two sessions.



MÖ PlasmaFill Spec:

Power supply	AC220V 50Hz
voltage:	110-220V, 50HZ
Max Temp:	90°C
Max Capacity	10 Syringes
Display	LED Digital
Power	200W
Machine size	290*280*140mm
Weight	4 KG
Packing weight	9.5KG

SCANDINAVIAN
SKINCARE SYSTEMS

Please contact us for price enquires, demonstrations and any other information:

info@scandinavianskincare.com
www.scandinavianskincare.com

Norway +47 951 00 442
UK +44 (0) 7711 086227